



### **HEINE 2.5x/420mm HR Binocular Loupe with 3S LED HeadLight®: Kit 2**

- 3S LED HeadLight® with headband battery (rechargeable) - mPack 30,000 Lumens at 250mm with 3SLED LoupeLight.
- 100% mobility - for absolute freedom of movement! Plug-in transformer. Direct operation possible while charging through the use of the mobile plug-in transformer. Charging time 2 hours. Top safety standards. Integrated charge and performance checks during operation. UL-compliant. System compatibility with mPack for even greater capacity.
- HR Technology for an even bigger field of view. The Brightest Image. Our exclusive coating technology reduces light reflection to 0.3 %, providing the brightest image of any comparable loupe. The Crispest Image with Color Correction. State-of-the-art HEINE optics with achromatic lenses provide a standard for image clarity, crispness and color compared with conventional Galilean loupes. Super-Wide, Super-Deep Field of View. Improvements in our HR (High Resolution) optical technology provide a super-wide, fully-useable, distortion-free field of view (up to 250 mm) and the greatest depth of field (up to 290 mm) to eliminate eye strain even over extended periods of use.
- Dust- and Water-Resistant. All HR loupes are dust and water-resistant to IP65/IEC 60529 (2000) to ensure ease of cleaning/disinfection, and years of trouble-free use.
- Flexibility for Precise Adjustment. HEINE Binocular loupes have independently adjustable optics left and right to select the exact inter-pupillary distance (PD). Conventional designs using a single adjuster are always a compromise. HR loupes can be adjusted precisely for each user for optimum optical performance by means of 5 different mounting positions on either frame or headband (with angled adaptor).
- Super Light. HR optics weigh only 42 grams for user-comfort even over extended periods. 2x Magnification: Ideal for ENT, Audiology, dental applications and general surgery, the 2x loupes provide our largest field of view and the best depth of field to ensure the image is in focus over an exceptionally-wide range. 2.5x Magnification: The increase in magnification over 2x is greater than the figures suggest. HR 2.5x loupes offer an excellent balance between high magnification, field of view, depth of field, and ease of use. The 2.5x loupes are ideal for a wide range of applications such as ENT, Audiology, dental, ophthalmic, cosmetic and reconstructive surgery.
- The HR 2.5x is our most popular model. Mounted on Professional L headband with S-Guard. Corrective frame can be purchased separately and your Rx can be clipped onto the S-Guard. Kit has plug in transformer and one mPack battery.