



INTELLIGENT HEARING

S Y S T E M S



The Video Visual Reinforcement Manual Audiometry System. It is Portable, Compact and Provides Unlimited Reinforcement Options.

Video VRA



Easy-to-Use Controls

Activate video reinforcers at the touch of a button on the video response box or the number keys on your computer keyboard. Easily change between left and right monitors by clicking on the space bar; or set the response box to control both right and left monitors using the response box buttons. Quickly change the video sets by scrolling through different video play-lists that you design using the keyboard up and down arrow buttons.

Fast and Simple Setup and Operation

Simply assign a video clip to a response box button by selecting a video file from a list of available video clips. you can copy and use your own video files to make more reinforcer choices available. Design your own play-lists for fast reinforcer switching at test time. It's simple and quick to change the reinforcing time period, sound volume. and to assign video clips to play on a right side or a left side monitor.



Unlimited Reinforcer Variety

Choose from the animations that are provided with your system, or use any MPEG or AVI movie of your preference. Use sets of clips that will appeal to the particular population the test subject belongs to. Create multiple video sets, and have them ready to go, to keep the patients full attention during testing.

VVRA Specifications

- Connects to any Windows® computer with video play capabilities and a free USB port.
- Multiple monitor capability in computers equipped with additional monitor outputs.
- Eight available video windows
 - Use four as right and four as left with dual monitors
 - Use all eight windows together with a single monitor



REV53467